

23 Gauge Two Step Trocar System

Reusable Two Step Trocar System with closure valves, 23Ga



5173-1

Package includes:

- Trocar Cannula with closure valves - 3 pcs
- Universal Infusion Line - 1 pc
- Blunt Cannula Inserter - 3 pcs
- Fixation Plate - 1 pc
- Loading Forceps - 1 pc



Reusable Two Step Trocar System with closure plugs, 23Ga

5173

Package includes:


- Trocar Cannula - 3 pcs
- Trocar Cannula Plug - 3 pcs
- Universal Infusion Line - 1 pc
- Blunt Cannula Inserter - 3 pcs
- Fixation Plate - 1 pc
- Loading Forceps - 1 pc



Watch video clip on [MicroVision website](#)

5186

Loading Forceps




5187

Blunt Cannula Inserter



5188

Fixation Plate



5202

Universal Infusion Line

Packaged sterile, five per box



23 Gauge Vitrectomy Accessoreis



Disposable Trocar Cannula Set, 23Ga

5190

Each Set includes:

- Instrument Cannula - 2ea
- Cannula Plug - 2ea

Sterile, five sets per box



Closure Valves for 23Ga Trocar Cannula

5209

Replacement Closure Valves - 1 set of 4pcs, sterile




5184 Angled MVR Blade, 23Ga - 6 per box



Self-Retaining Sutureless Cannulas, 23Ga

- Sutureless fixation
- Retaining ridge holds cannula in place
- Stainless steel cannula hub




5144 Self-Retaining Sutureless Infusion Cannula 23Ga, 4.0mm

5221 Self-Retaining Sutureless Infusion Cannula 23Ga, 6.0mm

5165 Self-Retaining Sutureless Silicone Oil Cannula 23Ga, 4.0mm


5222 Self-Retaining Sutureless Silicone Oil Cannula 23Ga, 6.0mm



Dual Bore PFC Cannula, 23Ga

5203 Dual Bore PFC Cannula, 23Ga - 5 per box

- Simultaneous infusion of heavy liquids and aspiration of intraocular fluids
- Dual bore enables to control constant intraocular pressure during injection



23 Gauge Two Step Trocar System

Disposable Two Step Trocar System, 23Ga



Trocar Cannula Inserter works together with Trocar Knife to accomplish Self-Sealing Trocar Cannula placement without using fixation plate.



5181 Two Step Trocar System, 23Ga

Each package includes:


- Three Trocar Cannula Inserters with preloaded Self-sealing Trocar Cannulas
- One Trocar Knife
- One Infusion Line

Not Available until further notice

Watch video clip on MicroVision website

Trocar Cannula Inserter

Trocar Cannula Inserter comes with Cannula pre-installed. Delivered in retracted position.




Self-Sealing Trocar Cannula

Stainless steel Trocar Cannula is complete with self-sealing cap, eliminating need for plugs.



Trocar Knife

Trocar Knife consists of sharp trocar mounted to plastic handle. Delivered in retracted position.




23 Gauge Vitrectomy Cutters

CE

Posterior Vitrectomy Cutters, 23Ga

5172 30 PSI	Alcon: Accurus - 2500 cpm DORC: Associate - 2500 cpm
5216	DORC: Harmony Total TTC - 1500 cpm
5171 25 PSI	B&L: Millennium, Premiere - 750 cpm Carl Zeiss Meditec/IOL Tech: Pentasys - 1500cpm Fritz Ruck: Pentasys - 1500 cpm Croma Pharma/Corneal: Open Phaco - 1500 cpm Optikon: Pulsar ² , Assistant - 2000 cpm Nidek: CV-12000, CV-6000, CV-7000 - 800 cpm Nidek: CV-24000 - 2500 cpm



Anterior Vitrectomy Cutters, 23Ga

5208 20 PSI	Alcon: Legacy, Universal AMO: Sovereign, Signature, AMO Plus, Prestige Staar Surgical: Wave Nidek: CV-12000, CV-6000, CV-7000
-----------------------	--


Fiber Optic Probes, 23Ga

5177	Fiber Optic Probe, 23Ga - Alcon Accurus
5178	Wide Angle Fiber Optic Probe, 23Ga - Alcon Accurus
5179	Fiber Optic Probe, 23Ga - B&L
5180	Wide Angle Fiber Optic Probe, 23Ga - B&L
5163	Fiber Optic Probe, 23Ga - Ocutome 5mm connector
5169	Wide Angle Fiber Optic Probe, 23Ga - Ocutome 5mm connector
5210	Shielded Fiber Optic Probe, 23Ga - Ocutome 5mm connector Panoramic Illumination with a centered focused spot




Silicone Tip Cannulas, 23Ga


5162	Brush Tip Cannula, 23Ga - 0.6mm silicone brush tip
5161	Soft Tip Cannula, 23Ga - 2.0mm silicone tip
5192	Diamond Dusted Soft Tip Cannula, 23Ga - 1.0mm silicone tip (for triamcinolone removal)
5164	Charles Flute Cannula, 23Ga - Blunt tip



Brush Tip



Soft Tip



Charles Flute