


25 Gauge Vitrectomy Products



One Step Trocar System, 25Ga

Insert Handle consists of sharp trocar mounted to plastic handle. Trocar cannula preloaded on Insert Handle as shown. Trocar Cannula has a sealing silicone valve, eliminating need for plugs.

Insert Handle with preloaded Trocar Cannula



5150 Trocar System, 25Ga

Each package includes:

- Four insert handles with Trocar Cannulas as shown
- One Infusion Line

Infusion Line



Trocar Cannula




Watch video clip on [MicroVision website](#)

Vitrectomy Cutters, 25Ga

5157
30 PSI **Alcon:** Accurus - 2500 cpm
DORC: Associate - 2500 cpm

B&L: Millennium, Premiere - 750 cpm
DORC: Harmony Total TTC - 1500 cpm
5022 - Pressure Line Adapter required

5140
25 PSI **Carl Zeiss Meditec/IOL Tech:** Pentasys - 1500 cpm
Fritz Ruck: Pentasys - 1500 cpm
Croma Pharma/Corneal: Open Phaco - 1500 cpm
Optikon: Pulsar², Assistant - 2000 cpm
AMO: Gemini - 2000 cpm
Nidek: CV-24000 - 2500 cpm
Nidek: CV-12000, VT-5000 - 800 cpm



5022 - Pressure Line Adapter - 5 per pack
Connects 5140 with DORC Harmony Total TTC

Fiber Optic Probes, 25Ga

- 5135** Fiber-Optic Probe, 25Ga - Alcon Accurus
5138 Fiber-Optic Probe, 25Ga - B&L
5126 Fiber-Optic Probe, 25Ga - Ocutome connector



Silicone Tip Cannulas, 25Ga

- 5160** Brush Tip Cannula, 25Ga - 0.6mm silicone brush tip
5152 Soft Tip Cannula, 25Ga - 2.0mm silicone tip
5156 Charles Flute Cannula, 25Ga - Blunt tip

